

# MC 780 PORTABLE

## SEGMENTAL MULTI FREQUENCY BODY COMPOSITION ANALYSER

### 3 FREQUENCY - SEGMENTAL BODY COMPOSITION REPORT IN SECONDS

#### MEASUREMENTS

- Body fat %
- Fat mass kg
- Fat free mass kg
- Muscle mass kg
- Total body water %
- Extra cellular water kg
- Intra cellular water kg
- Body mass index
- Bone mass kg
- Physique rating
- Visceral fat rating
- ECW/TBW ratio
- Basal metabolic rate kcal
- Phase angle
- Basal metabolic rate indicator
- Metabolic age
- Segmental readings for each leg, arm and trunk:  
Fat percentage, Fat distribution rating, Muscle mass rating, Muscle mass kg, Muscle mass balance, Leg muscle score, Body fat distribution and Reactance/resistance

#### SPECIFICATIONS

- Uses non-invasive 8 electrode multi-frequency Bio-electrical Impedance Analysis
- Print an in depth assessment report to any PictBrigde printer
- Weight capacity 270kg (max) x 0.1kg
- USB Connectivity
- Results can be easily saved to an SD Card
- Lightweight, easy to transport and set-up for use
- MDD class IIa and NAWI class III compliant

#### MATERIALS CODES

MC780MS01A	TEU MC-780MAS KG UK NL
MC780MS02A	TEU MC-780MAS KG DE FR



**NEW**